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**berkeleyhealth**  
[support@berkeleyhealth.com](mailto:support@berkeleyhealth.com)  
[www.berkeleyhealth.com](http://www.berkeleyhealth.com)

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
Do not open the sealed foil pack until you are ready to perform the test. Take out the test card and put it on a flat surface.

**Using the lancet:** Do not use the lancet if its extension has been removed or damaged. Wash your hands well with soap and warm water (in order to increase blood flow to your finger) and dry thoroughly.

- The lancet consists of a housing, containing the needle. The trigger which looks like a lever is on the upper side of the lancet and releases the needle. The front end of the lancet is hidden and protected by a protective cap. Remove the protective cap by performing one complete rotation of the cap. During the rotation, do not press down on the lever, otherwise the needle will be released. (Fig. 1)
- Once the cap is loose remove it and put aside. The lancet is ready for use now. (Fig. 2)
- Press the opening of the lancet against your middle fingertip and depress the lever until you hear a click. (Fig. 3)

**Applying a large drop of blood on the test field:** In just 3 minutes you will get to know about your cholesterol level. Please follow the instructions.

- Now, with your whole arm hanging down and your pricked finger pointing downwards, wait patiently until a big drop of blood has formed. It may be necessary to squeeze on both sides of your finger as illustrated. (Fig. 4)
- Allow one big drop of blood to fall on the test field. Ensure that your fingertip itself does not touch the test field. If insufficient blood is used or the test field is touched, the separation membrane might not work properly. (Fig. 5)
- At least as much blood as illustrated here has to be applied to the test field. **Now wait 3 minutes.** Test readings may be too high if you wait too long, or too low if you do not wait long enough. (Fig. 6)

**While waiting you can cover the pricked finger with the provided plaster:** Remove the protective foil from the plaster without touching the absorbent pad. Please take care that both skin areas for the adhesive parts of the plaster are dry, clean, fat free and uninjured and that the injured skin area is completely covered by the absorbent pad.

**Reading and interpretation of test results:**

- Then gently pull the cover off the card starting where indicated with the black bar. The focus point will be displayed as soon as the cover is completely removed. (Fig. 7)
- Now within 30 seconds move the colour chart until it shows the colour or colour intensity that corresponds best to the colour of the focus point. (Fig. 8) If you find it difficult to determine the matching colour field, apply the principle of elimination determining which colour fields do not match, until only one field is left. **The result will be less reliable if you need more than 30 seconds to match the colour of the focus point to one of the colour fields.**
- Turn the whole card over to find your estimated cholesterol level displayed in mg/dL or mmol/L in the left respectively right window in the guides of the colour chart. Which unit is better for me? The test result can be read off in mg/dL (milligram per deciliter) and in mmol/L (millimol per litre). You have a free choice of unit. To be able to compare future results, you should note down the value obtained and the corresponding unit. (Fig. 9)

**Normal readings:** For healthy adults, your cholesterol level is most probably in the **NORMAL** range if your reading on the colour chart shows values up to 5.2 mmol/l (200 mg/dl) or up to 3.9 mmol/l (150 mg/dl) for those at high risk.

**High levels:** Your cholesterol level is most probably **HIGH** if your reading shows values above 5.2 mmol/l (200 mg/dl) for healthy adults, or 3.9 mmol/l (150 mg/dl) for those at high risk. This is also the case if you have had difficulty deciding between the 5.2 mmol/l and the 5.9 mmol/l colour range. If level is elevated, visit your doctor soon for a more accurate determination of your cholesterol levels.

**The result indicator is not uniform in colour. There are darker and lighter green areas. How do I interpret the test?**  
Choose the areas which, taken together, are larger. **NOTE:** Small reddish-brownish areas indicate that erythrocytes (red blood cells) came into contact with the result indicator. Ignore these areas when interpreting the test result.

A reagent pad that is white or cream coloured indicates an **INVALID** test result.

Possible reasons for an invalid test result include a damaged foil pack, improper storage or a mistake when performing the test or a cholesterol level at or below 3.2 mmol/l (125 mg/dl). Please keep all test components and contact the berkeleyhealth.

**Important note: Always consult your doctor before making any important medical decisions!**

If your cholesterol level is higher than the indicated normal range: \*

You should take the result seriously, although a one-off raised level is not particularly meaningful. If your cholesterol level appears high, you should see your doctor, who will repeat the test using another method which may involve distinguishing between HDL (high-density cholesterol) and LDL (low-density cholesterol) as well as the triglyceride levels in your blood, in order to better assess your risk of arteriosclerosis. It might be helpful to take these instructions for use with you to your doctor to give a better idea of the test you have performed yourself.

\* As a general guideline, total cholesterol levels in healthy adults should be 5.0 mmol/l or less, according to UK National Health Service (2022).

**Product information:** This test uses the enzymes cholesterol esterase, cholesterol oxidase and peroxidase to activate a dye resulting in a colour change proportional to the cholesterol level in your blood.

The test components should be considered potentially hazardous materials, but represent no danger if you use all test components according to the instructions for use.

- In-vitro diagnostic medical device for self-testing
- Keep out of reach of children
- For external use only
- Do not use after expiry date
- Do not use a damaged lancet
- Do not dismantle the test card
- Do not open or ingest the enclosed desiccant
- Dispose of the test components and the desiccant in your normal household waste
- Store at 4-30°C

**Performance characteristics: \*\***

Clinical trials were conducted using this test to provide evidence of performance. Reference cholesterol tests used were a Reflotron whole blood dry chemistry analyser or a standard clinical laboratory cholesterol assay on venous blood. The following results were obtained from 227 pooled samples.

\*\* Brent A. Blumenstein, PhD, 1991 (Study not publicly available, please contact the manufacturer)

**Diagnostic Sensitivity: 92,4 %**

Percentage of true positives out of all individuals with a condition or disease by defining a positive test as a result of 5.2 mmol/l (200 mg/dl) or greater

**Diagnostic specificity: 81,3 %**

Percentage of true negatives out of all individuals who do not have a disease or condition by defining a positive test as a result of 5.2 mmol/l (200 mg/dl) or greater

**Linearity: 93,8 %** (interpretation by technician), **90,7 %** (self-interpreted) Agreement with reference cholesterol testing

**Positive predictive value: 87,1 %**

Percentage of true positives out of all positive tests by defining a positive test as a result of 5.2 mmol/l (200 mg/dl) or greater

**Negative predictive value: 88,6 %**

Percentage of true negatives out of all negative tests by defining a positive test as a result of 5.2 mmol/l (200 mg/dl) or greater

**Accuracy: 87,7 %**

Percentage of agreed results out of all tests

**Interference factors - the test result will be falsified by:**  
- Haemoglobin > 20 g/dl and bilirubin > 4.2 mg/dl  
- In excessive amounts, ascorbic acid (vitamin C), alpha-methylldopa or Novalgin lead to apparently low values  
- The haematocrit value influences the result in reverse dependence  
- Steroids such as epandrosterone, dehydroepiandrosterone, camposterol or sitosterol falsify the result if present in unusual concentrations.

**When to perform the test:** Ideally your cholesterol level should be measured on an empty stomach. This means that you should not have had anything to eat or drink for a period of 12 hours before carrying out the test. The test results should also be interpreted in daylight or under a good light source in order to be able to better assess the colour variations.

**Test components:** Instructions for use, foil pack with test card and desiccant, lancet, plaster.  
**Preparations:** You also need a clean paper tissue and a clock with a second hand or a timer.

**CE 0483** This product fulfils the requirements of the DIRECTIVE 79/98/EC OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 27 October 1998 on in vitro diagnostic medical devices and is classified as IVD for self-testing. A conformity assessment procedure was carried out by the notified body mdc (0483).

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